

LabVIEW Installation instructions

* Log on as an administrator or as a user with administrator privileges.

Installation -

=====

1. Place the LabVIEW installation DVD in your DVD Drive and wait for the opening screen, or invoke the installation, then click 'Install NI SOFTWARE'
2. Wait for the Installer to initialize and then click 'Next'
Click the radio button: "I only want to evaluate products", then click 'Next'.
3. Select the relevant products from the Evaluation Product List, then click "Next".
4. At the prompt click 'Yes' (recommended and optional) and "Next"
5. In the 'Destination Directory' screen accept the default or choose a location where you have enough disk space available. Click "Next"
6. License Agreement screen: Read and Accept and click 'Next' twice.
7. Choose your preferred option (Yes/No) for the "NI Update Service"
8. Wait for the installation to complete.
9. Click 'Finish' -> ' Restart Computer '

Network License setup -

=====

In All Programs

National Instruments -> NI license Manager

on the left top, click the arrow

choose - Manage Volume License Servers

enter the following license server address:

labview.technion.ac.il:27000

Click 'OK'